


## Firesafe™ - making oxygen therapy safer for everyone

4	EN	Firesafe™ Cannula Valve - Instructions for Use
5	DE	Firesafe™ Rückschlagventil – Bedienungsanleitung
6	FR	Vanne de canule Firesafe™ - Mode d'emploi
7	ES	Válvula para cánulas Firesafe™ - Instrucciones de uso
8	IT	Valvola per cannula Firesafe™ - Istruzioni per l'uso
9	NL	Firesafe™ Cannula Valve - Gebruiksaanwijzing
10	PT	Válvula de cânula Firesafe™ - Instruções de utilização
11	SV	Firesafe™ Slangventil - Användarinstruktioner
12	CZ	Kanylový ventil Firesafe™ - Pokyny k použití
13	PL	Zawór kaniuli Firesafe™ - Instrukcja obsługi
14	RU	Клапан для каниули Firesafe™ - Инструкция по эксплуатации
15	RO	Valva pentru canulă Firesafe™ - Instrucțiuni de utilizare
16	FI	Firesafe™ kanyyliventtiili - Käyttöohje
17	DA	Firesafe™ Cylinderventil - Brugsvejledning
18	SK	Kanylový ventil Firesafe™ - Návod na použitie
19	NO	Firesafe™ kanyleventil - Bruksanvisning
20	SL	Ventil Firesafe™ za kanilo - Navodila za uporabo
21	LT	Firesafe™ kaniulės vožtuvas - Naudojimo instrukcija
22	LV	Firesafe™ kanulas vārsts - Lietošanas instrukcija
23	HU	Használati útmutató Firesafe™ - kaniülselepe
24	EL	Βαλβίδα Κάνουλας Firesafe™ - Οδηγίες Χρήσης
25	ET	Firesafe™ kanüüli klapp - Kasutusjuhend
26	BG	Канюла с клапан Firesafe™ - Инструкции за употреба
27	HR	Ventil kanile Firesafe™ - Upute za uporabu
28	TR	Firesafe™ Kanül Valfi - Kullanım Talimatları

- EN Figure 1
- DE Abbildung 1
- FR Figure 1
- ES Figura 1
- IT Figura 1
- NL Afbeelding 1
- PT Figura 1
- SV Figur 1
- CZ Obrázku 1
- PL Rysunku 1
- RU Рисунок 1
- RO Figura 1
- FI Kuvan 1
- DA Figur 1
- SK Obrázku č. 1
- NO Figur 1
- SL Sliki 1
- LT Pav 1
- LV Attēlā 1
- HU Ábra 1
- EL Σχήμα 1
- ET Joonisel 1
- BG Фигура 1
- HR Slici 1
- TR Şekil 1

## Intended Use

The Firesafe™ Cannula Valve is a thermal fuse designed to extinguish an oxygen delivery tube fire and stop the flow of oxygen if the tube is accidentally ignited.

It is recommended that at least two Firesafe™ devices are fitted to each oxygen delivery circuit. The first should be positioned close to the source of supply and the second close to the patient as shown in figure 1. The Firesafe™ Cannula Valve provides no protection if a fire burns through the oxygen delivery tubing upstream of the device.

## Warnings

1. Read through this instruction before installing a Firesafe™ Cannula Valve. As with all medical equipment, attempting to use or install this device without a thorough understanding of its operation and limitations may result in patient or user injury.
2. This device is intended for use in oxygen delivery tubing that has a low static operating pressure not exceeding that stated in the device specification. It must not be used for other applications.
3. The device must be oriented correctly in the patient circuit as shown in figure 1. Failure to fit the device in the correct orientation will render it ineffective in the event of a fire.
4. This device offers a resistance to flow; consider the impact of this impedance on the system when installing the device.
5. Do not store or install this device near an open flame or near a source of excessive heat that is likely to exceed that stated in the Device Specification.
6. Keep this device free from oil and/or grease. In particular, do not lubricate the hose barb connectors.
7. Never administer oxygen or undertake oxygen therapy while smoking or when near an open flame.
8. Additional risk control measures (e.g. SpO2 monitoring with an alarm) are essential for patients who might suffer injury or death in the event that their oxygen flow is stopped unexpectedly.
9. The device is intended for single patient use and should not be moved between installations due to the risk of cross contamination.
10. This device is designed for use in the delivery of oxygen or oxygen enriched air. Do not use the device with any other gas.

## Cleaning

Clean the exterior surfaces of the device using an alcohol or disinfectant wipe. Do not immerse the device in any fluid or allow fluid to enter through the hose barb connections.

## Maintenance & Disposal

The Firesafe™ Cannula Valve is maintenance free and has an intended life of 4 years. Dispose of any Firesafe™ Cannula Valve 4 years after the manufacturing date stated on the device label. Once actuated, the Firesafe™ Cannula Valve cannot be reset and must be discarded.

## Device Specification

Resistance to Flow	≤ 1 kPa @ 2 l/min
	≤ 2 kPa @ 5 l/min
	≤ 16 kPa @ 15 l/min
Flow Rate <sup>(1)</sup>	0.7 to 20 l/min
Maximum Operating Pressure <sup>(2)</sup>	450 kPa
Maximum internal leakage following activation	10 ml/min
Maximum external leakage following activation	5 ml/min
Operating Temperature	0 to 50 °C (32 to 122 °F)
Transit & Storage Temperature range	-20 to 60 °C (-4 to 140 °F)
Humidity Range	0 to 100% RH
Gas type	Oxygen/oxygen enriched air

<sup>(1)</sup> The Firesafe™ Cannula Valve can be safely used at lower flow rates but may not meet the internal leakage rate specified when activated in the event of a fire

<sup>(2)</sup> Maximum static pressure delivered by the gas supply source

## Installation Instructions

**Warning!** When actuated, the Firesafe™ Cannula Valve will stop the gas flow. The pressure upstream of the device will rise until it reaches the maximum output pressure of the supply system. In order to ensure that the connection between the Firesafe™ Cannula Valve and the supply system can withstand this pressure, it is recommended that a type test be carried out with each style of trumpet connector or cannula tubing used.

**Note:** Due to the differences in the flexibility of the materials used for tubing and trumpet connectors and variations of inner bore diameter, it may not always be possible to fully engage the third barb as shown in figure 1.

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BPR Medical Limited